

NEW PRODUCTS - AAO, 2021



TID Pharos Illuminated Depressor


- TID Pharos is our second generation transscleral illuminator developed with Dr. Bozho Todorich
- The complete, one-piece illuminated depressor connects directly into your vitrectomy machine
- Conventional ball-point depressor tip with a “Side Fire Illumination” feature allowing for a better view, less glare and improved surgeon ergonomics
- The TID Pharos is designed to offer the surgeon the ability to simultaneously perform unassisted scleral depression while shaving the vitreous base

ORDERING INFORMATION:

VS0801A - TID Pharos Illuminated Depressor - Alcon Constellation

VS0801B - TID Pharos Illuminated Depressor - Bausch and Lomb Stellaris

TID PHAROS ILLUMINATED DEPRESSOR



30ga PC Maintainer Supplementary Infusion Line



- Hands-free infusion line designed to help prevent hypotony while suturing
- Perfect for re-ops and repeat vitrectomies
- 30ga One-Step insertion needle with stability footplate for easy insertion
- Ideal for fellowship programs and teaching facilities
- Developed with Dr. Justin Gutman

ORDERING INFORMATION:

VS0240.30 - 30ga P.C. Maintainer Supplementary Infusion Line

Myopic ACTU8 Forceps

- Additional 5mm length added to the forceps (30mm vs 35mm) without compromising the benefits and advantages of the ACTU8 Forceps Product Line
- Not offered in a disposable version by Alcon or DORC

ORDERING INFORMATION:

VS0742.23-35 — 23ga Tewari (T2) Myopic Forceps

VS0742.25-35 — 25ga Tewari (T2) Myopic Forceps

VS0744.23-35 — 23ga ILM Myopic Forceps

VS0744.25-35 — 25ga ILM Myopic Forceps

ACTU8
premium disposable instrumentation



Vortex Surgical
680 Crown Industrial Ct., Suite F
Chesterfield, MO 63005

NEW PRODUCTS - AAO, 2021




27ga MAXReach Laser Probe

- 27ga addition to our popular MAXReach laser probe line
- 90+ degree curve for greater peripheral access
- Maximized for energy delivery on the DORC EVA
- Direct connect to all major US laser systems

ORDERING INFORMATION:

VS0140.27 - 27ga MAXReach Laser Probe



End-Grip ACTU8 Forceps



- Tapered grasping platforms provide utility of use for different applications
- Tweezer style forceps design
- Textured platforms to increase holding strength

ORDERING INFORMATION:

VS0743.23 - 23ga ACTU8 End-Grip Forceps

VS0743.25 - 25ga ACTU8 End-Grip Forceps

Blunt Tip Cannulas

- Radiused/Rounded tip to allow for use on or around the retinal surface
- Increase flow rates - 23ga - 5X, 25ga - 3X and 27ga - 20X

ORDERING INFORMATION:

VS0201.23 - 23ga Blunt Tip Cannula

VS0201.25 - 25ga Blunt Tip Cannula

VS0201.27 - 27ga Blunt Tip Cannula




27ga Side Port Cannulas

- Designed for delivery of staining agents, providing the surgeon more control over delivery of the stain
- Permits for a controlled injection that covers the surface of the membrane
- Available in single and dual port configurations


ORDERING INFORMATION:

VS0202.27 - 27ga Side Port Cannula

VS0203.27 - 27ga Dual Port Cannula



Side Port



Dual Port