


EZpass

SOFT TIP CANNULA



- Soft Tip Cannulas to fit any procedure
- Custom tubing for high flow performance
- Easy passage through all cannula systems without cutting the tip

1mm and 0.5mm
Soft Tip lengths available

BLUNT CANNULAS

- Highly radiused distal tip makes safe for retinal surface
- Flow rate increased significantly over soft tips
- Available in 23ga, 25ga and 27ga



SIDE PORT CANNULAS

- Unique closed end, side port opening allowing for controlled injections of stains and other therapies
- Available in 27ga single and dual side ports



SUBRETINAL INJECTION

- 41ga x 88mm polyimide tubing at tip
- Rigid polyimide allows for easy penetration through pathology
- Perfect for TPA injections
- Available in 25ga



VERSATILE SOFT TIP

- Extendable silicone tubing allows for atraumatic subretinal drainage.
- Extendable tubing length up to 11.5mm
- Available in 23ga and 25ga



BACKFLUSH CANNULA

- High Flow Soft Tip with backflush reflux handle
- Maximized volume reflux button
- Available in 23ga, 25ga and 27ga
- Retractable and Non-Retractable versions available



Catalog Numbering and Ordering information are located on the reverse side





Catalog Numbers and Descriptions

SOFT TIP CANNULA

VS0200.23	23ga Soft Tip - 1mm, Sterile, 5/bx
VS0200.23B	23ga Soft Tip - 1mm, Sterile, 100/bx
VS0200.25	25ga Soft Tip - 1mm, Sterile, 5/bx
VS0200.25B	25ga Soft Tip - 1mm, Sterile, 100/bx
VS0200.27	27ga Soft Tip - 1mm, Sterile, 5/bx
VS0205.23	23ga Soft Tip - 0.5mm, Sterile, 5/bx
VS0205.23B	23ga Soft Tip - 0.5mm, Sterile, 100/bx
VS0205.25	25ga Soft Tip - 0.5mm, Sterile, 5/bx
VS0205.25B	25ga Soft Tip - 0.5mm, Sterile, 100/bx
VS0205.27	27ga Soft Tip - 0.5mm, Sterile, 5/bx

BLUNT CANNULA

VS0201.23	23ga Blunt, Sterile, 5/bx
VS0201.25	25ga Blunt, Sterile, 5/bx
VS0201.27	27ga Blunt, Sterile, 5/bx

BACKFLUSH CANNULA

VS0270.23	23ga Backflush, Sterile, 5/bx
VS0270.25	25ga Backflush, Sterile, 5/bx
VS0270.27	27ga Backflush, Sterile, 5/bx
VS0275.23	23ga Ret. Backflush, Sterile, 5/bx
VS0275.25	25ga Ret. Backflush, Sterile, 5/bx
VS0275.27	27ga Ret. Backflush, Sterile, 5/bx

SIDE PORT CANNULA

VS0202.27	27ga Side Port-Single, Sterile, 100/bx
VS0203.27	27ga Side Port-Dual, Sterile, 100/bx

SUBRETINAL INJECTION CANNULA

VS0220.25	25ga Subretinal Injection, Sterile, 5/bx
-----------	--

VERSATILE SOFT TIP CANNULA

VS0280.23	23ga Versatile Soft Tip, Sterile, 5/bx
VS0280.25	25ga Versatile Soft Tip, Sterile, 5/bx

ORDERING INFORMATION:



CANNULAS

Vortex Surgical
4 Research Park Drive, Suite 124
Saint Charles, MO 63304 USA

Phone: 636-778-4350
Email: info@vortexsurgical.com
Web: www.vortexsurgical.com



VORTEX SURGICAL