

BEACON

ILLUMINATION SYSTEM

A complete illumination system designed for the Ocular Oncologist for use in the clinic or operating room!




BEACON LED LIGHT SOURCE VS0900

- 50 watt LED Light Source
- Customized illumination pathway to maximize illumination
- Small and Portable for use in the clinic or operating room

BEACON TRANSILLUMINATOR VS0802, sterile, 5/bx

- A conventional ball-point depressor tip with a “Side Fire Illumination” feature allows for a better view, less glare and improved surgeon ergonomics
- The device is designed to offer the surgeon the ability to simultaneously perform unassisted scleral depression while performing oncology evaluations
- Additional benefits include: identifying retained lens fragments, identifying radiation seed placements and locating breaks or tears in the outer periphery.



PATENT PENDING



www.vortexasurgical.com